

LCORL Polyclonal antibody

Catalog Number: 17468-1-AP

Basic Information

Catalog Number:

17468-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11377

GenBank Accession Number:

BC037322

GeneID (NCBI):

254251

UNIPROT ID:

Q8N3X6

Full Name:

ligand dependent nuclear receptor
corepressor-like

Calculated MW:

318 aa, 35 kDa, 67 kDa

Observed MW:

67 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2400

IHC 1:20-1:200

Applications

Tested Applications:

IHC, WB, ELISA

Species Specificity:

human, mouse, rat

Note-IHC: suggested antigen retrieval with
**TE buffer pH 9.0; (*) Alternatively, antigen
retrieval may be performed with citrate
buffer pH 6.0**

Positive Controls:

WB : mouse testis tissue, PC-3 cells

IHC : human lung tissue,

Background Information

Storage

Storage:


Store at -20°C. Stable for one year after shipment.

Storage Buffer:


PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 17468-1-AP (LCORL antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human lung using 17468-1-AP (LCORL antibody) at dilution of 1:100 (under 10x lens).