

For Research Use Only

ITGBL1 Polyclonal antibody

Catalog Number: 17484-1-AP

Featured Product

3 Publications



Basic Information

Catalog Number:

17484-1-AP

Concentration:

190 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11521

GenBank Accession Number:

BC036788

GeneID (NCBI):

9358

UNIPROT ID:

O95965

Full Name:

integrin, beta-like 1 (with EGF-like repeat domains)

Calculated MW:

494 aa, 54 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC, IF

Species Specificity:

human, mouse, rat

Cited Species:

human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Zhe Wang	35596183	Biomed Eng Online	WB, IHC
Qing Ji	32139701	Nat Commun	WB, IHC, IF
Li Sun	27261588	Biomed Pharmacother	WB, IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

