

For Research Use Only

PDS5A Polyclonal antibody, PBS Only

Catalog Number: 17485-1-PBS

Featured Product



Basic Information

Catalog Number:

17485-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11525

GenBank Accession Number:

BC009650

GeneID (NCBI):

23244

UNIPROT ID:

Q29RF7

Full Name:

PDS5, regulator of cohesion maintenance, homolog A (S. cerevisiae)

Calculated MW:

1297 aa, 147 kDa

Observed MW:

130 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 17485-1-AP (PDS5A antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 17485-1-PBS in a different storage buffer formulation.



WB result of PDS5A antibody (17485-1-AP; 1:4000; incubated at room temperature for 1.5 hours) with sh-Control and sh-PDS5A transfected HEK-293 cells. This data was developed using the same antibody clone with 17485-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human normal colon slide using 17485-1-AP (PDS5A antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 17485-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed A431 cells using PDS5A antibody (17485-1-AP) at dilution of 1:400 and CoraLite@488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-phalloidin (red). This data was developed using the same antibody clone with 17485-1-PBS in a different storage buffer formulation.