

For Research Use Only

# CFD Polyclonal antibody

Catalog Number: 17486-1-AP

1 Publications



## Basic Information

Catalog Number:

17486-1-AP

Size:

187 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11548

GenBank Accession Number:

BC057807

GeneID (NCBI):

1675

UNIPROT ID:

P00746

Full Name:

complement factor D (adipsin)

Calculated MW:

253 aa, 27 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

IF

Species Specificity:

human, mouse

Cited Species:

human

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Binshan Zhang	36582241	Front Immunol	IF

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

