

For Research Use Only

ICMT Polyclonal antibody

Catalog Number: 17511-1-AP



Basic Information

Catalog Number:

17511-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_012405

GeneID (NCBI):

23463

UNIPROT ID:

O60725

Full Name:

isoprenylcysteine carboxyl
methyltransferase

Calculated MW:

32 kDa

Observed MW:

32-34 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : Jurkat cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells lysates were subjected to SDS PAGE followed by western blot with 17511-1-AP (ICMT antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.