

For Research Use Only

CTF1-Specific Polyclonal antibody

Catalog Number: 17518-1-AP



Basic Information

Catalog Number:

17518-1-AP

Size:

150 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001330

GeneID (NCBI):

1489

UNIPROT ID:

Q16619

Full Name:

cardiotrophin 1

Calculated MW:

21 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human heart tissue slide using 17518-1-AP (CTF1-Specific Antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human heart tissue slide using 17518-1-AP (CTF1-Specific Antibody) at dilution of 1:50 (under 40x lens).