

For Research Use Only

# PLEKHA9 Polyclonal antibody

Catalog Number: 17568-1-AP



## Basic Information

Catalog Number:

17568-1-AP

Size:

800 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11772

GenBank Accession Number:

BC063575

GeneID (NCBI):

51054

UNIPROT ID:

O95397

Full Name:

pleckstrin homology domain  
containing, family A  
(phosphoinositide binding specific)  
member 9

Calculated MW:

391 aa, 44 kDa

Observed MW:

65 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : Jurkat cells, HeLa cells, MCF-7 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17568-1-AP (PLEKHA9 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.