

For Research Use Only

Inhibin Beta B Polyclonal antibody

Catalog Number: 17577-1-AP

Featured Product

7 Publications



Basic Information

Catalog Number:

17577-1-AP

Size:

240 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10956

GenBank Accession Number:

BC030029

GeneID (NCBI):

3625

UNIPROT ID:

P09529

Full Name:

inhibin, beta B

Calculated MW:

407 aa, 45 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Cai Kailai K	21998670	PLoS One	WB
Min Tao	31631336	J Fish Biol	IHC
Mohamed M'baye	26063610	J Reprod Dev	

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

