

For Research Use Only

# Inhibin Beta B Polyclonal antibody

Catalog Number: 17577-1-AP **7 Publications**



## Basic Information

<b>Catalog Number:</b> 17577-1-AP	<b>GenBank Accession Number:</b> BC030029	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 240 µg/ml	<b>GeneID (NCBI):</b> 3625	
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P09529	
<b>Isotype:</b> IgG	<b>Full Name:</b> inhibin, beta B	
<b>Immunogen Catalog Number:</b> AG10956	<b>Calculated MW:</b> 407 aa, 45 kDa	

## Applications

**Tested Applications:**  
ELISA

**Cited Applications:**  
WB, IHC

**Species Specificity:**  
human, mouse, rat

**Cited Species:**  
human, mouse

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Cai Kailai K	21998670	PLoS One	WB
Min Tao	31631336	J Fish Biol	IHC
Mohamed M'baye	26063610	J Reprod Dev	

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data