

For Research Use Only

ATP5H Polyclonal antibody

Catalog Number: 17589-1-AP

6 Publications



Basic Information

Catalog Number:

17589-1-AP

Concentration:

400 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11429

GenBank Accession Number:

BC032245

GeneID (NCBI):

10476

UNIPROT ID:

O75947

Full Name:

ATP synthase, H⁺ transporting, mitochondrial F₀ complex, subunit d

Calculated MW:

137 aa, 16 kDa

Observed MW:

19-22 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:12000

IP 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

IHC 1:50-1:500

IF/ICC 1:50-1:500

Applications

Tested Applications:

WB, IHC, IF/ICC, IP, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : Jurkat cells, mouse liver tissue, rat brain tissue, rat liver tissue

IP : mouse liver tissue,

IHC : human lung cancer tissue, human pancreas tissue

IF/ICC : HepG2 cells,

Background Information

Mitochondrial membrane ATP synthase (F₁-F₀ ATP synthase or Complex V) produces ATP from ADP in the presence of a proton gradient across the membrane which is generated by electron transport complexes of the respiratory chain. It is composed of the soluble catalytic core, F₁, and the membrane-spanning component and F₀, which comprises the proton channel. The F₀ seems to have nine subunits (a, b, c, d, e, f, g, F₆ and 8). ATP5H gene encodes ATP synthase subunit d of the F₀ complex.

Notable Publications

Author	Pubmed ID	Journal	Application
Liangjun Xia	35401830	Theranostics	WB
Xia Feng	35867854	Brain	WB
Margalida A Frau-Méndez	27338255	Brain Pathol	IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Jurkat cells were subjected to SDS PAGE followed by western blot with 17589-1-AP (ATP5H antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using 17589-1-AP (ATP5H antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



IP result of anti-ATP5H (IP:17589-1-AP, 3ug; Detection:17589-1-AP 1:2000) with mouse liver tissue lysate 6000ug.



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using ATP5H antibody (17589-1-AP) at dilution of 1:200 and CoraLite®594-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-4), CL488-Phalloidin (green).