

For Research Use Only

ASCC3 Polyclonal antibody, PBS Only

Catalog Number:17627-1-PBS



Basic Information

Catalog Number:

17627-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11848

GenBank Accession Number:

BC050681

GeneID (NCBI):

10973

UNIPROT ID:

Q8N3C0

Full Name:

activating signal cointegrator 1
complex subunit 3

Calculated MW:

111 aa, 13 kDa, 251 kDa

Observed MW:

251 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human

Background Information

ASCC3 (Activating Signal Cointegrator 1 Complex Subunit 3) is a multifunctional protein involved in various cellular processes, including transcriptional regulation, DNA repair, and ribosome quality control. It is a member of the ASCC (Activating Signal Cointegrator 1 Complex) family, which functions as a bridge between co-repressors and co-activators in transcriptional regulation.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 17627-1-AP (ASCC3 antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 17627-1-PBS in a different storage buffer formulation.

IP result of anti-ASCC3 (IP:17627-1-AP, 4ug; Detection:17627-1-AP 1:400) with HeLa cells lysate 1080ug. This data was developed using the same antibody clone with 17627-1-PBS in a different storage buffer formulation.