

For Research Use Only

PLA2G3 Polyclonal antibody, PBS Only

Catalog Number:17646-1-PBS



Basic Information

Catalog Number:

17646-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11845

GenBank Accession Number:

BC025316

GeneID (NCBI):

50487

UNIPROT ID:

Q9NZ20

Full Name:

phospholipase A2, group III

Calculated MW:

509 aa, 57 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data