

For Research Use Only

# ENPEP Polyclonal antibody

Catalog Number: 17655-1-AP

3 Publications



## Basic Information

Catalog Number:

17655-1-AP

Concentration:

207 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11903

GenBank Accession Number:

BC094770

GeneID (NCBI):

2028

UNIPROT ID:

Q07075

Full Name:

glutamyl aminopeptidase  
(aminopeptidase A)

Calculated MW:

957 aa, 109 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF/ICC 1:10-1:100

## Applications

Tested Applications:

IF/ICC, ELISA

Cited Applications:

IF

Species Specificity:

human

Cited Species:

mouse, rat

Positive Controls:

IF/ICC : HepG2 cells,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Taha A Jan	34260939	Cell Rep	IF
Tae Jin Lee	19743427	Small	FC
Liwen Luo	39985222	Adv Sci (Weinh)	

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of HepG2 cells, using ENPEP antibody 17655-1-AP at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).