

For Research Use Only

# ARHGAP30 Polyclonal antibody

Catalog Number: 17706-1-AP



## Basic Information

|                                      |  |   |
|--------------------------------------|--|---|
| Catalog Number:<br>17706-1-AP        | GenBank Accession Number:<br>BC053688          | Purification Method:<br>Antigen affinity purification |
| Size:<br>N/A                         | GeneID (NCBI):<br>257106                       |   |
| Source:<br>Rabbit                    | UNIPROT ID:<br>Q7Z6I6                          |   |
| Isotype:<br>IgG                      | Full Name:<br>Rho GTPase activating protein 30 |   |
| Immunogen Catalog Number:<br>AG11992 | Calculated MW:<br>890 aa, 95 kDa               |   |

## Applications

Tested Applications:  
ELISA

Species Specificity:  
human, mouse, rat

## Background Information

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data