

For Research Use Only

ULBP1 Polyclonal antibody

Catalog Number: 17715-1-AP **10 Publications**



Basic Information

Catalog Number: 17715-1-AP	GenBank Accession Number: BC035416	Purification Method: Antigen affinity purification
Size: 700 µg/ml	GeneID (NCBI): 80329	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: Q9BZM6	
Isotype: IgG	Full Name: UL16 binding protein 1	
Immunogen Catalog Number: AG12079	Calculated MW: 244 aa, 28 kDa	
	Observed MW: 28 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HEK-293 cells, Jurkat cells
Cited Applications: WB, IF, IHC	
Species Specificity: human	
Cited Species: human, mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Jiahui Zhang	30268986	Mol Immunol	WB,IF
Madison N Temples	33463353	ACS Biomater Sci Eng	IF
Charlotte Flüh	29356965	Histochem Cell Biol	IHC

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 17715-1-AP (ULBP1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.