

For Research Use Only

DCP2 Polyclonal antibody

Catalog Number: 17744-1-AP



Basic Information

Catalog Number:

17744-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12001

GenBank Accession Number:

BC064593

GeneID (NCBI):

167227

UNIPROT ID:

Q8IU60

Full Name:

DCP2 decapping enzyme homolog (S. cerevisiae)

Calculated MW:

385aa, 44 kDa; 420aa, 48 kDa

Observed MW:

48 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HepG2 cells, Caco-2 cells, T-47D cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17744-1-AP (DCP2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.