

For Research Use Only

PARP8 Polyclonal antibody, PBS Only

Catalog Number:17754-1-PBS



Basic Information

Catalog Number:

17754-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12111

GenBank Accession Number:

BC075801

GeneID (NCBI):

79668

UNIPROT ID:

Q8N3A8

Full Name:

poly (ADP-ribose) polymerase family, member 8

Calculated MW:

812 aa, 91 kDa

Observed MW:

91-96 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IP, Indirect ELISA

Species Specificity:

human

Background Information

Protein mono-ADP-ribosyltransferase PARP8, also known as ADP-ribosyltransferase diphtheria toxin-like 16 (ARTD16), is primarily localised on the nuclear envelope for the majority of the cell cycle but localises to centrosomes and spindle poles during mitosis.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-PARP8 (IP:17754-1-AP, 4ug; Detection:17754-1-AP 1:500) with THP-1 cells lysate 2480 ug. This data was developed using the same antibody clone with 17754-1-PBS in a different storage buffer formulation.