

For Research Use Only

SPIRE2 Polyclonal antibody

Catalog Number: 17757-1-AP

2 Publications



Basic Information

Catalog Number:

17757-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12117

GenBank Accession Number:

BC111030

GeneID (NCBI):

84501

UNIPROT ID:

Q8WWL2

Full Name:

spire homolog 2 (Drosophila)

Calculated MW:

666 aa, 74 kDa

Observed MW:

79 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IP, IF, IHC

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Positive Controls:

WB : mouse cerebellum tissue,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Lixia Hao	36122882	Neuroscience	WB, IF
Yoon Jung Hae JH	22775139	J Proteome Res	WB, IP, IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse cerebellum tissue were subjected to SDS PAGE followed by western blot with 17757-1-AP (SPIRE2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.