

For Research Use Only

CACNB4 Polyclonal antibody

Catalog Number: 17770-1-AP **1 Publications**



Basic Information

Catalog Number: 17770-1-AP	GenBank Accession Number: BC075049	Purification Method: Antigen affinity purification
Size: 500 µg/ml	GeneID (NCBI): 785	Recommended Dilutions: WB 1:1000-1:4000 IHC 1:50-1:500
Source: Rabbit	UNIPROT ID: O00305	
Isotype: IgG	Full Name: calcium channel, voltage-dependent, beta 4 subunit	
Immunogen Catalog Number: AG12178	Calculated MW: 520 aa, 58 kDa	
	Observed MW: 58 kDa	

Applications

Tested Applications: IHC, WB, ELISA	Positive Controls:
Cited Applications: WB	WB : mouse brain tissue, rat brain tissue
Species Specificity: human, mouse, rat	IHC : mouse brain tissue,
Cited Species: human	

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Cui Ma	33230455	Mol Ther Nucleic Acids	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17770-1-AP (CACNB4 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 17770-1-AP (CACNB4 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 17770-1-AP (CACNB4 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).