

For Research Use Only

ORC6 Polyclonal antibody

Catalog Number: 17784-1-AP

Featured Product

3 Publications



Basic Information

Catalog Number:

17784-1-AP

Concentration:

650 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11769

GenBank Accession Number:

BC063565

GeneID (NCBI):

23594

UNIPROT ID:

Q9Y5N6

Full Name:

origin recognition complex, subunit 6 like (yeast)

Calculated MW:

252 aa, 28 kDa

Observed MW:

28 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:12000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB : HeLa cells, mouse spleen tissue, Jurkat cells, MCF-7 cells, MDA-MB-231 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Rui Mao	33188157	Aging (Albany NY)	WB, IHC
Lifei Fan	28467790	Oncotarget	WB
Jinfeng Zhu	37901236	Front Immunol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17784-1-AP (ORC6 antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours.