

For Research Use Only

# ORC6 Polyclonal antibody, PBS Only

Catalog Number: 17784-1-PBS

Featured Product



## Basic Information

Catalog Number:

17784-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11769

GenBank Accession Number:

BC063565

GeneID (NCBI):

23594

UNIPROT ID:

Q9Y5N6

Full Name:

origin recognition complex, subunit 6 like (yeast)

Calculated MW:

252 aa, 28 kDa

Observed MW:

28 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17784-1-AP (ORC6 antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 17784-1-PBS in a different storage buffer formulation.