

For Research Use Only

SLC25A33 Polyclonal antibody

Catalog Number: 17794-1-AP **1 Publications**



Basic Information

| | | |
|---|--|--|
| Catalog Number: 17794-1-AP | GenBank Accession Number: BC073135 | Purification Method: Antigen affinity purification |
| Size: 350 µg/ml | GeneID (NCBI): 84275 | Recommended Dilutions: WB 1:500-1:1000 |
| Source: Rabbit | UNIPROT ID: Q9BSK2 | |
| Isotype: IgG | Full Name: solute carrier family 25, member 33 | |
| Immunogen Catalog Number: AG11926 | Calculated MW: 321 aa, 35 kDa | |
| | Observed MW: 35-40 kDa | |

Applications

| | |
|--|---|
| Tested Applications: WB, ELISA | Positive Controls: WB : U2Os cells, |
| Cited Applications: WB | |
| Species Specificity: human | |
| Cited Species: human | |

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|--------------|-----------|----------------------------|-------------|
| Xiao-Jing Gu | 36438699 | World J Gastrointest Oncol | WB |

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



200kDa
150kDa
100kDa
75kDa
50kDa
25kDa
17kDa

anti-PRC1
anti-SL

200kDa
150kDa
100kDa
75kDa
50kDa
25kDa
17kDa

anti-PRC1 from PTCs in HeLa

WB result of anti-SLC25A33 (17794-1-AP) by Dr. Geraldine Kelly, University College Cork, Ireland.