

For Research Use Only

SLC25A33 Polyclonal antibody

Catalog Number: 17794-1-AP

1 Publications



Basic Information

Catalog Number:

17794-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11926

GenBank Accession Number:

BC073135

GeneID (NCBI):

84275

UNIPROT ID:

Q9BSK2

Full Name:

solute carrier family 25, member 33

Calculated MW:

321 aa, 35 kDa

Observed MW:

35-40 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : U2Os cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xiao-Jing Gu	36438699	World J Gastrointest Oncol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



WB result of anti-SLC25A33 (17794-1-AP) by Dr. Geraldine Kelly, University College Cork, Ireland.