

For Research Use Only

FBXO3 Polyclonal antibody

Catalog Number: 17803-1-AP



Basic Information

Catalog Number:

17803-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12071

GenBank Accession Number:

BC046110

GeneID (NCBI):

26273

UNIPROT ID:

Q9UK99

Full Name:

F-box protein 3

Calculated MW:

471 aa, 55 kDa

Observed MW:

48-50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, HepG2 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 17803-1-AP (FBXO3 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.