

For Research Use Only

RHOQ/TC10 Polyclonal antibody

Catalog Number: 17805-1-AP

Featured Product

4 Publications



Basic Information

Catalog Number:

17805-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12078

GenBank Accession Number:

BC070485

GeneID (NCBI):

23433

UNIPROT ID:

P17081

Full Name:

ras homolog gene family, member Q

Calculated MW:

205 aa, 23 kDa

Observed MW:

23 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue, human brain tissue, rat brain tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Shingo Koinuma	29072354	Genes Cells	
Shingo Koinuma	33156548	J Neurochem	WB
Marianne Burbage	28747344	J Immunol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17805-1-AP (RHOQ/TC10 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.