

For Research Use Only

# ARHGAP12 Polyclonal antibody

Catalog Number: 17855-1-AP



## Basic Information

|                                             |                                                       |                                                              |
|---------------------------------------------|-------------------------------------------------------|--------------------------------------------------------------|
| <b>Catalog Number:</b><br>17855-1-AP        | <b>GenBank Accession Number:</b><br>BC094719          | <b>Purification Method:</b><br>Antigen affinity purification |
| <b>Size:</b><br>N/A                         | <b>GeneID (NCBI):</b><br>94134                        |                                                              |
| <b>Source:</b><br>Rabbit                    | <b>UNIPROT ID:</b><br>Q8IWW6                          |                                                              |
| <b>Isotype:</b><br>IgG                      | <b>Full Name:</b><br>Rho GTPase activating protein 12 |                                                              |
| <b>Immunogen Catalog Number:</b><br>AG12309 | <b>Calculated MW:</b><br>841 aa, 96 kDa               |                                                              |

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human, mouse, rat

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data