

For Research Use Only

# GH1 Polyclonal antibody

Catalog Number: 17867-1-AP **2 Publications**



## Basic Information

Catalog Number:

17867-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12350

GenBank Accession Number:

BC075012

GeneID (NCBI):

2688

UNIPROT ID:

P01241

Full Name:

GH1

Calculated MW:

217 aa, 25 kDa

Observed MW:

22 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Hong Yao	31176743	Cancer Lett	WB
Lian Chaoqun C	24306150	Curr Microbiol	

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data