

For Research Use Only

MMRN1 Polyclonal antibody

Catalog Number: 17878-1-AP **1 Publications**



Basic Information

Catalog Number: 17878-1-AP	GenBank Accession Number: BC063848	Purification Method: Antigen affinity purification
Size: 400 µg/ml	GeneID (NCBI): 22915	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q13201	
Isotype: IgG	Full Name: multimerin 1	
Immunogen Catalog Number: AG11739	Calculated MW: 1228 aa, 138 kDa	
	Observed MW: 138-150 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Human peripheral blood platelets,
Cited Applications: WB	
Species Specificity: human	
Cited Species: mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Na Wu	36037246	Proteomics	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Human peripheral blood platelets were subjected to SDS PAGE followed by western blot with 17878-1-AP (MMRN1 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.