

For Research Use Only

MMRN1 Polyclonal antibody

Catalog Number: 17878-1-AP

1 Publications



Basic Information

Catalog Number:

17878-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11739

GenBank Accession Number:

BC063848

GeneID (NCBI):

22915

UNIPROT ID:

Q13201

Full Name:

multimerin 1

Calculated MW:

1228 aa, 138 kDa

Observed MW:

138-150 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

mouse

Positive Controls:

WB : Human peripheral blood platelets,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Na Wu	36037246	Proteomics	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Human peripheral blood platelets were subjected to SDS PAGE followed by western blot with 17878-1-AP (MMRN1 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.