

For Research Use Only

RALY Polyclonal antibody

Catalog Number: 17896-1-AP



Basic Information

Catalog Number:

17896-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12250

GenBank Accession Number:

BC103753

GeneID (NCBI):

22913

UNIPROT ID:

Q9UKM9

Full Name:

RNA binding protein, autoantigenic
(hnRNP-associated with lethal yellow
homolog (mouse))

Calculated MW:

307 aa, 33 kDa

Observed MW:

34-36 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells, HeLa cells, Jurkat cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17896-1-AP (RALY antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.