

For Research Use Only

Opticin Polyclonal antibody

Catalog Number: 17897-1-AP



Basic Information

Catalog Number:

17897-1-AP

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12257

GenBank Accession Number:

BC074943

GeneID (NCBI):

26254

UNIPROT ID:

Q9UBM4

Full Name:

opticin

Calculated MW:

37 kDa

Observed MW:

45 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue, mouse cerebellum tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 17897-1-AP (Opticin antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.