

For Research Use Only

LCN1 Polyclonal antibody

Catalog Number: 17900-1-AP

1 Publications



Basic Information

Catalog Number:

17900-1-AP

Size:

900 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12270

GenBank Accession Number:

BC074925

GeneID (NCBI):

3933

UNIPROT ID:

P31025

Full Name:

lipocalin 1 (tear prealbumin)

Calculated MW:

19 kDa

Observed MW:

17-19 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : human saliva tissue, Transfected HEK-293 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
John Mello Camille C Guzman	32977557	Biosensors (Basel)	

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human saliva were subjected to SDS PAGE followed by western blot with 17900-1-AP (LCN1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.