

For Research Use Only

NEUROG3 Polyclonal antibody

Catalog Number: 17902-1-AP

1 Publications



Basic Information

Catalog Number:

17902-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12272

GenBank Accession Number:

BC074938

GeneID (NCBI):

50674

UNIPROT ID:

Q9Y4Z2

Full Name:

neurogenin 3

Calculated MW:

23 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|--------------------|-----------|--------------|-------------|
| Priyadharsini Raju | 25661254 | Reproduction | WB |

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

