

For Research Use Only

AP1G2 Polyclonal antibody

Catalog Number: 17907-1-AP



Basic Information

Catalog Number:

17907-1-AP

Concentration:

550 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12300

GenBank Accession Number:

BC051833

GeneID (NCBI):

8906

UNIPROT ID:

O75843

Full Name:

adaptor-related protein complex 1,
gamma 2 subunit

Calculated MW:

640 aa, 71 kDa

Observed MW:

87 kDa, 67 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : THP-1 cells, mouse brain tissue, MCF-7 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



THP-1 cells were subjected to SDS PAGE followed by western blot with 17907-1-AP (AP1G2 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.