

For Research Use Only

# SOX4 Polyclonal antibody

Catalog Number: 17919-1-AP

Featured Product

3 Publications



## Basic Information

Catalog Number:

17919-1-AP

Size:

370 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12333

GenBank Accession Number:

BC072668

GeneID (NCBI):

6659

UNIPROT ID:

Q06945

Full Name:

SRY (sex determining region Y)-box 4

Calculated MW:

474 aa, 47 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Zhen-Zhen Chen	27553854	Nat Commun	WB, IHC
Rongfei Han	27446439	Oncol Lett	WB
Maisa S Abduh	38402910	Biochem Pharmacol	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

