

For Research Use Only

KCNN1 Polyclonal antibody

Catalog Number: 17929-1-AP



Basic Information

Catalog Number:

17929-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12357

GenBank Accession Number:

BC075037

GeneID (NCBI):

3780

UNIPROT ID:

Q92952

Full Name:

potassium intermediate/small conductance calcium-activated channel, subfamily N, member 1

Calculated MW:

543 aa, 60 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB: mouse brain tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 17929-1-AP (KCNN1 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.