

For Research Use Only

MIP Polyclonal antibody

Catalog Number: 17977-1-AP



Basic Information

Catalog Number:

17977-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12369

GenBank Accession Number:

BC074913

GeneID (NCBI):

4284

UNIPROT ID:

P30301

Full Name:

major intrinsic protein of lens fiber

Calculated MW:

263 aa, 28 kDa

Observed MW:

26-28 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse eye tissue, mouse skeletal muscle tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse eye tissue were subjected to SDS PAGE followed by western blot with 17977-1-AP (MIP antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.