

For Research Use Only

ADH1A Polyclonal antibody

Catalog Number: 17998-1-AP **1 Publications**



Basic Information

Catalog Number: 17998-1-AP	GenBank Accession Number: BC074738	Purification Method: Antigen affinity purification
Size: N/A	GeneID (NCBI): 124	
Source: Rabbit	UNIPROT ID: P07327	
Isotype: IgG	Full Name: alcohol dehydrogenase 1A (class I), alpha polypeptide	
Immunogen Catalog Number: AG12467	Calculated MW: 375 aa, 40 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
IHC

Species Specificity:
human, mouse, rat

Cited Species:
human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Karim Sumera S	23300006	Toxicol Sci	IHC

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data