

For Research Use Only

# GREM1/Gremlin 1 Polyclonal antibody, PBS Only

Catalog Number: 18024-1-PBS



## Basic Information

Catalog Number:

18024-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12588

GenBank Accession Number:

BC093778

GeneID (NCBI):

26585

UNIPROT ID:

O60565

Full Name:

gremlin 1, cysteine knot superfamily,  
homolog (Xenopus laevis)

Calculated MW:

184 aa, 21 kDa

Observed MW:

21 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 18024-1-AP (GREM1/Gremlin 1 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 18024-1-PBS in a different storage buffer formulation.