

For Research Use Only

# KCNA1 Polyclonal antibody

Catalog Number: 18067-1-AP



## Basic Information

Catalog Number:

18067-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12699

GenBank Accession Number:

BC101733

GeneID (NCBI):

3736

UNIPROT ID:

Q09470

Full Name:

potassium voltage-gated channel,  
shaker-related subfamily, member 1  
(episodic ataxia with myokymia)

Calculated MW:

495 aa, 56 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

IF/ICC, ELISA

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data