

For Research Use Only

# SIGLEC10 Polyclonal antibody

Catalog Number: 18091-1-AP



## Basic Information

Catalog Number:

18091-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12814

GenBank Accession Number:

BC101725

GeneID (NCBI):

89790

UNIPROT ID:

Q96LC7

Full Name:

sialic acid binding Ig-like lectin 10

Calculated MW:

697 aa, 77 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

