

For Research Use Only

OTOR Polyclonal antibody, PBS Only

Catalog Number: 18132-1-PBS



Basic Information

Catalog Number:

18132-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12714

GenBank Accession Number:

BC101688

GeneID (NCBI):

56914

UNIPROT ID:

Q9NRC9

Full Name:

otoraplin

Calculated MW:

128 aa, 14 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human liver cancer tissue slide using 18132-1-AP (OTOR antibody) at dilution of 1:100 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 18132-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissue slide using 18132-1-AP (OTOR antibody) at dilution of 1:100 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 18132-1-PBS in a different storage buffer formulation.