

For Research Use Only

# TNFSF9 Polyclonal antibody

Catalog Number: 18140-1-AP



## Basic Information

Catalog Number:

18140-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12768

GenBank Accession Number:

BC104805

GeneID (NCBI):

8744

UNIPROT ID:

P41273

Full Name:

tumor necrosis factor (ligand)  
superfamily, member 9

Calculated MW:

254 aa, 27 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

TNFSF9, also named as 4-1BBL, belongs to the tumor necrosis factor family. It is a cytokine that binds to TNFRSF9. TNFSF9 induces the proliferation of activated peripheral blood T-cells. It may have a role in activation-induced cell death (AICD) and may play a role in cognate interactions between T-cells and B-cells/macrophages.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

