

For Research Use Only

RAB9B Polyclonal antibody

Catalog Number: 18152-1-AP



Basic Information

Catalog Number:

18152-1-AP

Size:

200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12848

GenBank Accession Number:

BC093756

GeneID (NCBI):

51209

UNIPROT ID:

Q9NP90

Full Name:

RAB9B, member RAS oncogene family

Calculated MW:

201 aa, 23 kDa

Observed MW:

23-25 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

IF/ICC, WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, mouse brain tissue, HepG2 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 18152-1-AP (RAB9B antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.



IF result of RAB9B antibody (18152-1-AP, 1:200) detected with PFA fixed HeLa cells. Green: RAB9B antibody stains late endosomes. Red, Golgi stain.