

For Research Use Only

# OMP Polyclonal antibody

Catalog Number: 18185-1-AP **1 Publications**



## Basic Information

<b>Catalog Number:</b> 18185-1-AP	<b>GenBank Accession Number:</b> BC069365	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 350 µg/ml	<b>GeneID (NCBI):</b> 4975	<b>Recommended Dilutions:</b> WB 1:200-1:1000
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P47874	
<b>Isotype:</b> IgG	<b>Full Name:</b> olfactory marker protein	
<b>Immunogen Catalog Number:</b> AG12899	<b>Calculated MW:</b> 163 aa, 19 kDa	
	<b>Observed MW:</b> 20 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : mouse brain tissue,
<b>Cited Applications:</b> WB	
<b>Species Specificity:</b> human, mouse	

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Jian Wu	34297096	Int Immunol	WB


## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 18185-1-AP (OMP Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.