

For Research Use Only

DYSFIP1 Polyclonal antibody, PBS Only

Catalog Number: 18208-1-PBS



Basic Information

| | | |
|---|--|--|
| Catalog Number: 18208-1-PBS | GenBank Accession Number: BC093707 | Purification Method: Antigen affinity purification |
| Concentration: 1 mg/ml | GeneID (NCBI): 116729 | |
| Source: Rabbit | UNIPROT ID: Q86WC6 | |
| Isotype: IgG | Full Name: dysferlin interacting protein 1 | |
| Immunogen Catalog Number: AG12839 | Calculated MW: 154 aa, 17 kDa | |
| | Observed MW: 17 kDa | |

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human, rat

Background Information

DYSFIP1 is also known as PPP1R27. DYSFIP1 inhibits phosphatase activity of protein phosphatase 1 (PP1) complexes. The molecular weight of DYSFIP1 is 17 kDa.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Rat skeletal muscle tissue were subjected to SDS PAGE followed by western blot with 18208-1-AP (DYSFIP1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 18208-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed THP-1 cells using DYSFIP1 antibody (18208-1-AP) at dilution of 1:200 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 18208-1-PBS in a different storage buffer formulation.