

For Research Use Only

HIST1H4F Polyclonal antibody

Catalog Number: 18341-1-AP



Basic Information

Catalog Number:

18341-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13113

GenBank Accession Number:

BC093763

GeneID (NCBI):

8361

UNIPROT ID:

P62805

Full Name:

histone cluster 1, H4f

Calculated MW:

103 aa, 11 kDa

Observed MW:

11 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : human spleen tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human spleen tissue were subjected to SDS PAGE followed by western blot with 18341-1-AP (HIST1H4F antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.