

For Research Use Only

# SMAD6 Polyclonal antibody

Catalog Number: 18342-1-AP **1 Publications**



## Basic Information

Catalog Number: 18342-1-AP	GenBank Accession Number: BC012986	Purification Method: Antigen affinity purification
Size: N/A	GeneID (NCBI): 4091	
Source: Rabbit	UNIPROT ID: O43541	
Isotype: IgG	Full Name: SMAD family member 6	
Immunogen Catalog Number: AG13128	Calculated MW: 53 kDa	

## Applications

Tested Applications:  
ELISA

Cited Applications:  
WB

Species Specificity:  
human, mouse, rat

Cited Species:  
human

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Wei Lu	33174391	J Cell Mol Med	WB

## Storage

Storage:  
Store at -20°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data