

For Research Use Only

IGLV6-57 Polyclonal antibody

Catalog Number: 18344-1-AP



Basic Information

Catalog Number:

18344-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13141

GenBank Accession Number:

BC023973

GeneID (NCBI):

28778

Full Name:

immunoglobulin lambda variable 6-57

Calculated MW:

148 aa, 16 kDa

Observed MW:

25 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human plasma,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Human plasma was subjected to SDS PAGE followed by western blot with 18344-1-AP (IGLV6-57 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.