

For Research Use Only

SPATA20 Polyclonal antibody

Catalog Number: 18373-1-AP

1 Publications



Basic Information

Catalog Number:

18373-1-AP

Size:

200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13234

GenBank Accession Number:

BC025255

GeneID (NCBI):

64847

UNIPROT ID:

Q8TB22

Full Name:

spermatogenesis associated 20

Calculated MW:

88 kDa

Observed MW:

80-88 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF, CoIP

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : rat testis tissue,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xiang Wang	36415156	Clin Genet	WB, IF, CoIP

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat testis tissue were subjected to SDS PAGE followed by western blot with 18373-1-AP (SPATA20 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.